



CERTIFICATION by CUSTODIAN of RECORDS

Tucson Police Department Crime Laboratory
270 S. Stone Avenue Tucson, AZ 85701-1917

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number 80-002265 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

TS4
9.19.2011

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.

Date: 08/25/11 Cal 1 0.099 AC Cal 2 0.097 AC

Between the dates of 07/28/11 and 08/25/11 the instrument was operating properly and accurately.

The instrument was taken out of service for repair on N/A

The instrument was returned to service on N/A

Cheryl K. Straub

Cheryl K. Straub


09/08/2011

ADPS QAS Permit # 33408

SUBSCRIBED AND SWORN to before me this 8 day of Sept, 20 11.

Ann L. Brecroft

Notary Public



*Includes documents specified in A.A.C. R13-10-104.A.5

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

187 CS

EXHIBIT G-4

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

TS
9.19.2011

QA SPECIALIST Cheryl Straub AGENCY Tucson Police Department
DATE 08/25/2011 TIME 14:05
INTOXILYZER SERIAL # 80-002265 LOCATION CRIME LAB

(✓) 1. Display reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

(✓) 1. Clock time check.
(✓) 2. Date check.

OPERATIONAL TESTS

(✓) 1. Alcohol-free subject test result 0. 000 AC.
(✓) 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
(✓) 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
(✓) 4. Standard Calibration Check standard 0. 100 AC
Result: 0. 097 AC

Instrument is operating properly and accurately. YES ✓ NO

COMMENTS

SIGNATURE

Cheryl Straub

287
cg

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number : 80-002265

08/25/2011

14:07:00

134

Standard Lot#: 20910100A1

Last Changed By: C. STRAUB

#50772

QAS: C. STRAUB

#50772

TUCSON PD

Operator: C. STRAUB

#50772

TUCSON PD

Subject: SQAP, SQAP

DOB: 11/11/1111

Sex: F

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:07:55
Diagnostic Test	Pass	14:08:23
Air Blank	0.000	14:08:50
0.100 Cal Check	0.099	14:09:09
Air Blank	0.000	14:09:38
Subject Test	0.000	14:10:42
Air Blank	ABT*	14:10:50
Air Blank	0.000	14:11:18

*Sequence Aborted

Not a Successfully Completed Test Sequence

377
CS

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number : 80-002265

08/25/2011

14:11:22

BSH

Standard Lot#: 20910100A1

Last Changed By: C. STRAUB

#50772

QAS: C. STRAUB

#50772

TUCSON PD

Operator: C. STRAUB

#50772

TUCSON PD

Subject: SQAP, SQAP

DOB: 11/11/1111

Sex: F

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	14:11:59
Air Blank	IPS**	14:12:07

*RFI Detect

**Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number : 80-002265

08/25/2011

14:12:13

Standard Lot#: 20910100A1

Last Changed By: C. STRAUB

#50772

QAS: C. STRAUB

#50772

TUCSON PD

Operator: C. STRAUB

#50772

TUCSON PD

Subject: SQAP, SQAP

DOB: 11/11/1111

Sex: F

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:13:05
Diagnostic Test	Pass	14:13:32
Air Blank	0.000	14:14:00
0.100 Cal Check	0.099	14:14:18
Air Blank	0.000	14:14:48
Subject Test	0.019*	14:17:56
Air Blank	0.000	14:18:31
Air Blank	0.000	14:18:58
Subject Test	XXX**	14:19:26
Air Blank	0.000	14:19:59

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

577
CS

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number: 80-002265

08/25/2011

14:20:19

TS4


Standard Lot#: 20910100A1

Last Changed By: C. STRAUB #50772

QAS: C. STRAUB #50772

TUCSON PD

Test	g/210L	Time
Air Blank	0.000	14:20:46
Cal Check	0.097	14:21:05
Air Blank	0.000	14:21:35



Operator's Signature

697
08

EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

B4
9.19.2011

QA Specialist Cheryl STRAUB AGENCY Tucson Police Department
DATE 08/25/2011 TIME 14:05
Intoxilyzer Serial # 80-002265 Location CRIME LAB

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. 100 AC.
OR
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 099 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Cheryl Straub

DPS Form Exh G-2 (Rev. 05-01)

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number: 80-002265

08/25/2011

14:05:11

787
05
Bh


Standard Lot#: 20910100A1

Last Changed By: C. STRAUB #50772

QAS: C. STRAUB #50772

TUCSON PD

Test	g/210L	Time
Air Blank	0.000	14:05:38
Cal Check	0.099	14:05:57
Air Blank	0.000	14:06:27



Operator's Signature